

+ ALLEVYN[®] Ag NON-ADHESIVE Antimicrobial Foam Dressing

Smith+Nephew

Helping you get **CLOSER TO ZERO[®]**
human and economic consequence of wounds



+ ALLEVYN[◇] Ag NON-ADHESIVE Foam Dressing with advanced triple layer construction and addition of silver sulfadiazine to the foam pad to protect the wound from invasive pathogenic microorganisms.^{1,3}

Features and benefits

Helps minimise maceration

ALLEVYN Ag NON-ADHESIVE Dressing maintains a moist wound environment when exudate is present to help minimise the risk of maceration.⁴⁻⁷

Bacterial barrier

ALLEVYN Ag NON-ADHESIVE Dressing provides an effective barrier to microbial contamination to protect the wound from invasive pathogenic microorganisms shown in-vitro.¹⁻³

Easy to use

ALLEVYN Ag NON-ADHESIVE Dressing is easy to use, apply and remove and can be cut to size.^{4,8-10}

Suitable for use under compression

ALLEVYN Ag NON-ADHESIVE Dressing may be used under compression bandages and retains fluid when compressed.^{5,12,13}

Indications

ALLEVYN Ag NON-ADHESIVE Hydrocellular Foam Dressing is an antimicrobial absorbent dressing for the management, by secondary intention of chronic and acute full thickness, partial thickness or shallow granulating, exuding wounds such as:

- Pressure injuries
- Venous ulcers
- Diabetic ulcers
- Burns
- Donor sites
- Fungating /malignant wounds
- Surgically dehisced wounds

ALLEVYN Ag NON-ADHESIVE may be used on infected wounds. Where the product is used on infected wounds the wounds should be treated as per local protocol.

For detailed product information, including indications for use, contraindications, precautions and warnings, please consult the product's applicable Instructions for Use (IFU) prior to use.

ALLEVYN Ag NON-ADHESIVE order references

| Sizes | Code | Carton |
|-------------|----------|--------|
| 5cm x 5cm | 66800083 | 10 |
| 10cm x 10cm | 66800086 | 10 |
| 15cm x 15cm | 66800089 | 10 |
| 20cm x 20cm | 66800091 | 10 |

Smith & Nephew Pty Ltd
Australia
T +61 2 9857 3999
F +61 2 9857 3900

smith-nephew.com/australia

Smith & Nephew Ltd
New Zealand
T +64 9 820 2840
F +64 9 820 2841

smith-nephew.com/new-zealand

◇Trademark of Smith+Nephew
All Trademarks acknowledged
© 2020 Smith+Nephew
25338 V1-anz 0920

References

1. Smith+Nephew 2007. Antimicrobial activity of ALLEVYN Ag NON-ADHESIVE Dressing against a broad spectrum of microorganisms. Internal report. DOF 0703006.
2. Smith+Nephew 2007. Bacterial barrier testing (Wet-Wet) on samples of ALLEVYN Ag NON-ADHESIVE Dressing with a 7 day test duration against *S. marcescens*. Internal report. DOF 070304.
3. Smith+Nephew 2007. Bacterial barrier testing (Wet-Wet) on samples of ALLEVYN Ag NON-ADHESIVE Dressing with a 7 day test duration against MRSA. Internal report. DOF 070305.
4. Smith & Nephew 2006. A prospective user evaluation on the effect of process changes to the manufacture of ALLEVYN NON-ADHESIVE wound dressings. Internal report CER/CE/022/ANA.
5. Smith+Nephew 2009. ALLEVYN Ag NON-ADHESIVE Dressing physical properties. Internal report. DS/09/013/R8.
6. Smith+Nephew 2008. A multi-centre in-market evaluation of ALLEVYN Ag Dressings. Internal report. SR/CIME/009.
7. Smith+Nephew 2019. Use of Moisture Vapour Pemeability (MVP) and Moisture Vapour Transmission Rate (MVTR) data to support product claims referring to moist wound healing. Internal statement. EO.AWM.PCSgen.001.v2.
8. Foster AVM, Greenhill MT, Edmonds ME. Comparing two dressings in the treatment of diabetic foot ulcers. *Journal of Wound Care*. 1994;3(5):224-228.
9. Smith+Nephew 2015. Cutting of ALLEVYN variants. Internal report DS/14/318/R.
10. Kurring PA, Roberts CD, Quinlan D. Evaluation of a hydrocellular dressing in the management of exuding wounds in the community. *British Journal of Nursing* 1994;3(20):1049-1050, 1052-1043.
11. Leonard S, McCluskey P, Long S, et al. An evaluation of ALLEVYN ADHESIVE and NON-ADHESIVE Foam Dressings. *Wounds UK*. 2009;5(1):17-28.
12. Smith+Nephew 2018. Project Peridot testing - ALLEVYN NON-ADHESIVE Dressings. Internal report. DS/18/264/R.
13. Franks PJ, Moody M, Moffatt CJ, et al. Randomized trial of two foam dressings in the management of chronic venous ulceration. *Wound Repair Regen*. 2007;15(2):197-202.